



---

# **Development of an Automated Homogenizer –Autogizer and Its Application in Brain Uptake Studies**

**Shiyong Wang<sup>1</sup>, Hong Mei<sup>1</sup>, Kwokei Ng<sup>1</sup>, Keith Workowski<sup>2</sup>,  
Tom Astle<sup>2</sup> and Walter Korfmacher<sup>1</sup>,**

**Department of Drug Metabolism and Pharmacokinetics, <sup>1</sup>Schering-Plough  
Research Institute, 2015 Galloping Hill Road, Kenilworth, NJ 07033 and  
<sup>2</sup>Tomtec, Inc. Hamden, CT 06514**



# INTRODUCTION

---

## BACKGROUND

- The evaluation of brain uptake of drug candidates is essential in drug discovery.
- Sample preparation for brain tissue sample is labor intensive and limits the turnover time of the analysis.
- Long-term exposure to noise from manual homogenization of tissue samples is a health hazard.

## OBJECTIVE

- To design an automated homogenizer to process brain samples at much faster rate
- To improve the working environment




# Picture of Autogizer





# Method: Brain Sample Homogenization



15-mL centrifuge tube

Homogenize  
use Autogizer




## Parameters setting

- OSC speed: 4
- cycle: 3
- cycle time: 45 sec.
- 3 washes. each wash 35 sec.

**Total time: 2.5 min/5 samples**



# Design for Washing Reservoir




**Direction of water  
flow**

**Using continuously up-  
flowing water to rinse the  
probe is much more  
efficient than rinsing the  
tube in static water**




# Experiment Design: Carry-Over of Autogizer



## Washing reservoir:

Tubes in Column A, B, C were filled with continuously flowing water from bottom



## Sample Processing Tray:

Column 1: Blank brain tissue


Column 2: Brain tissue spiked with 5000 ng/g of 4 Standards

Column 3: Blank brain tissue

Column 4: Blank brain tissue

## Procedures:


Autogizer was programmed to move with the following order:



Aliquots from 1, 3, and 4 were taken for LC/MS analysis.



# Brain Homogenate Preparation for LC/MS Analysis






# Result: Carry-over of Autogizer after Processing 5000 ng/g Brain Standards

	Concentration (ng/g)			
	CMPD_A	CMPD_B	CMPD_C	CMPD_D
<b>Column 1 (before processing 5000 ng/g standards)</b>				
tube 1	0	0	0	0
tube 2	0	1	0	0
tube 3	0	1	0	0
tube 4	0	1	1	0
tube 5	0	0	0	0
<b>Column 3 (after processing 5000 ng/g standards and 3 washes with water)</b>				
<b>tube 1</b>	<b>1</b>	<b>0</b>	<b>0</b>	<b>3</b>
<b>tube 2</b>	<b>0</b>	<b>0</b>	<b>0</b>	<b>1</b>
<b>tube 3</b>	<b>1</b>	<b>0</b>	<b>1</b>	<b>3</b>
<b>tube 4</b>	<b>1</b>	<b>0</b>	<b>0</b>	<b>1</b>
<b>tube 5</b>	<b>1</b>	<b>0</b>	<b>0</b>	<b>1</b>
<b>Column 4 (after Column 3)</b>				
tube 1	0	0	0	1
tube 2	0	0	0	0
tube 3	0	0	0	1
tube 4	0	0	0	1
tube 5	0	0	0	1



## Results: Calibration Curve of Brain S standards



	Front Standards		Back Standards	
Conc.ng/g	Back-calc.	%DEV	Back-calc.	%DEV
1	1	20	0.9	-9.7
2.5	2	-2.2	2	-12
5	6	9.2	6	24
10	10	-0.6	11	12
25	24	-3.9	23	-9.2
50	60	20	56	13
100	106	5.8	95	-5.5
250	257	2.9	266	6.5
500	600	20	514	2.9
1000	875	-13	989	-1.1
2500	2457	-1.7	2367	-5.3
5000	4794	-4.1	5538	11





## Results: Within-run Reproducibility in Brain Tissue

	<b>Conc. (ng/g)</b>	<b>Calculated Conc. (ng/g)</b>	<b>Mean (N=6) (ng/g)</b>	<b>CV %</b>	<b>Accuracy</b>
<b>1</b>	<b>5</b>	<b>5.5; 5.9; 4.6; 4.6; 6.1; 5</b>	<b>5</b>	<b>12</b>	<b>106%</b>
<b>2</b>	<b>50</b>	<b>53; 55; 60; 49; 38; 64</b>	<b>53</b>	<b>18</b>	<b>106%</b>
<b>3</b>	<b>500</b>	<b>453; 436; 531; 499; 468; 524</b>	<b>485</b>	<b>8</b>	<b>97%</b>



# Results: An Example of Brain Uptake Study





## Results: Comparison of Manual Homogenization and Auto Homogenization (based on 20 samples)

	Manual	Auto
<b>Total Time</b>	100 min	10 min
<b>Labor Time</b>	100 min for a set of brains	~0. No labor needed. The robot did the job
<b>Noise Hazard</b>	Yes. Person exposed to noise for 40 min	No. The robot was set in a separated room.
<b>Safety Compliance</b>	Person used his/her hand holding the vibrating probe for 40 min, prone to accidents	The robot did the job without human intervention
<b>Solvent Consumption</b>	20L of water and 2L of MEOH	4L of water



# CONCLUSIONS

---

- The designed Autogizer can process brain samples at least 10 times faster than manual homogenization.
- The precision and accuracy are very good and meet our acceptance criteria for drug discovery sample analysis
- It has  $< 0.1\%$  carry-over and consumes no organic solvent and uses 5 times less water to rinse the probe.
- Because it is automated and can be left alone, it is noise hazard free and more safety compliant.