

## Harvester Specifications

	<b>Model Numbers</b>	<b>Length</b>	<b>Width</b>	<b>Height</b>	<b>Weight</b>
<b>Mach II Manual</b>	96-IIM	20.5"	14.0"	19.0"	60 lbs.
<b>Mach II Automatic</b>	96-II	24.0"	14.0"	19.0"	60 lbs.
<b>Mach III Automatic</b>	96-III, IIIU, IIIW	24.0"	12.5"	15.5"	50 lbs.
<b>Mach III Manual</b>	96-IIIM, IIIMU, IIIMW	19.5"	12.5"	14.0"	50 lbs.
<b>Harvester 9600</b>	All Models	16.5"	16.5"	25.5"	80 lbs.
<b>Harvex</b>		19.5"	12.5"	14.0"	50 lbs.

Weight and dimensions are approximate.

## System Requirements

	<b>Electrical</b>	<b>Air Input Min/Max</b>	<b>Vacuum</b>
<b>All Manual Models</b>	100-125Vac, 50/60 Hz @ 1.0 Amp 220-240Vac, 50/60Hz @ 0.5 Amp	5/20 PSI	25 + in Hg @ 8 + CFM
<b>All Automatic Models</b>	100-125Vac, 50/60 Hz @ 1.0 Amp 220-240Vac, 50/60 Hz @ 0.5 Amp	5/20 PSI	25 + in Hg @ 8 + CFM
<b>Harvester 9600 All Models</b>	100-120Vac, 50/60 Hz @ 3.0 Amp 230-240Vac, 50/60 Hz @ 1.5 Amp	35/45 PSI	25 + in Hg @ 8 + CFM
<b>Harvex</b>	100-125Vac, 50/60 Hz @ 1.0 Amp 220-240Vac, 50/60Hz @ 0.5 Amp	5/20 PSI	25 + in Hg @ 8 + CFM