



Larger View



LIVE BOTTOM RECIRCULATING RESERVOIR

Use with cell cultures, magnetic particles and SPA beads. Specifically designed to keep particulate matter in a uniform suspension for aspiration by the pipettor. A peristaltic pump is included to keep the entire 3" x 5" reservoir in constant motion. The reservoir is made of Teflon to permit autoclaving and is inert to DMSO and most other chemicals commonly used in HTS.

Catalog No.: 196-370-31 (for use with the Quadra 3)

Catalog No.: 196-370- 51) (for use with the Quadra-Plus")

TOMTEC